

Anti-Phospho-EIF4B-Ser406 antibody (STJ196297)

STJ196297

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Eukaryotic translation initiation factor 4B-Ser406 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1975
Gene Symbol	EIF4B
Uniprot ID	IF4B_HUMAN
Immunogen	Synthesized phospho peptide around human eIF4B (Ser406)
Immunogen Region	
Specificity	This antibody detects endogenous levels of Human Mouse Rat eIF4B (phospho-Ser406)
Immunogen Sequence	