

Anti-FAT3 antibody (896-946 aa) (STJ196598)

STJ196598

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protocadherin Fat 3 (896-946 aa) is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	IHC-P 1:50-200 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	120114
Gene Symbol	FAT3
Uniprot ID	FAT3_HUMAN
Immunogen	Synthesized peptide derived from the human FAT3 at the amino acid range 896-946
Immunogen Region	896-946 aa
Specificity	This antibody detects endogenous levels of FAT3 at Human/Mouse/Rat
Immunogen Sequence	