

Anti-SARS-nCoV-2 antibody [S06FC] (STJ196635)

STJ196635

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Human monoclonal antibody anti-SARS-CoV-2 Nucleoprotein is suitable for use in ELISA, ECL and Immunogold research applications.
Applications	ELISA/ECL/Immunogold
Host/Source	Human
Reactivity	COVID-19

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	S06FC
Concentration	2 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from 293T expression by affinity-chromatography using specific immunogen.
Dilution Range	ELISA 1:5000-1:100000 WB 1:1000-1:10000
Formulation	Liquid in PBS
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	The recombinant SARS-CoV-2 (2019-nCoV) N Protein
Immunogen Region	
Specificity	This antibody detects SARS-CoV-2 (2019-nCoV) Nucleoprotein
Immunogen Sequence	