

## Anti-JAK2 antibody [2F6] (STJ196685)

STJ196685

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Tyrosine-protein kinase JAK2 is suitable for use in Immunohistochemistry research applications.
<b>Applications</b>	IHC
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	2F6
<b>Concentration</b>	1 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
<b>Dilution Range</b>	IHC 1:50-200
<b>Formulation</b>	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	3717
<b>Gene Symbol</b>	JAK2
<b>Uniprot ID</b>	JAK2_HUMAN
<b>Immunogen</b>	Synthesized peptide derived from the human JAK2
<b>Immunogen Region</b>	
<b>Specificity</b>	This antibody detects endogenous levels of JAK2 at Human, Mouse, Rat
<b>Immunogen Sequence</b>	