

## Anti-MS4A1 antibody [ABT219R] (STJ197055)

STJ197055

### GENERAL INFORMATION

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-B-lymphocyte antigen CD20 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	IHC/IF
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	ABT219R
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	The antibody was affinity-purified from Cell supernatant by affinity-chromatography using specific immunogen.
<b>Dilution Range</b>	IHC-P 1:10 IF 1:50-200
<b>Formulation</b>	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
<b>Isotype</b>	IgG
<b>Storage</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	<a href="#">931</a>
<b>Gene Symbol</b>	<a href="#">MS4A1</a>
<b>Uniprot ID</b>	<a href="#">CD20_HUMAN</a>
<b>Immunogen</b>	Synthesized peptide derived from the human CD20
<b>Immunogen Region</b>	
<b>Specificity</b>	This antibody detects endogenous levels of CD20 at Human
<b>Immunogen Sequence</b>	