

Anti-PTGS2 antibody [ABT31R] (STJ197059)

STJ197059

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Prostaglandin G/H synthase 2 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	ABT31R
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from Cell supernatant by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:100 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5743
Gene Symbol	PTGS2
Uniprot ID	PGH2_HUMAN
Immunogen	Synthesized peptide derived from the human COX-2
Immunogen Region	
Specificity	This antibody detects endogenous levels of COX-2 at Human
Immunogen Sequence	