

Anti-S100A6 antibody [ABT179R] (STJ197072)

STJ197072

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein S100-A6 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	ABT179R
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from Cell supernatant by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:2000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6277
Gene Symbol	S100A6
Uniprot ID	S10A6_HUMAN
Immunogen	Synthesized peptide derived from the human S100A6
Immunogen Region	
Specificity	This antibody detects endogenous levels of S100A6 at Human
Immunogen Sequence	