

Anti-POU2AF1 antibody (150-256 aa) [ABT121] (STJ197102)

STJ197102

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-POU domain class 2-associating factor 1 (150-256 aa) is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	ABT121
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC-P 1:100-500 IF 1:100-500
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1k
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5450
Gene Symbol	POU2AF1
Uniprot ID	OBF1_HUMAN
Immunogen	Synthesized peptide derived from the human BOB.1 at the amino acid range 150-256
Immunogen Region	150-256 aa
Specificity	The antibody can specifically recognize human BOB.1 protein.
Immunogen Sequence	