

Anti-Tri-Methyl-Histone H2B-Lys5 antibody (STJ197297)

STJ197297

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tri-Methyl-Histone H2B-Lys5 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS pH7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3017/8339/8343/8344/8346/8347 255626 H2BC4.H2BC6.H2BC7.H2BC8.H2BC10 H2BC1
Uniprot ID	H2B1C_HUMAN H2B1A_HUMAN H2B1B_HUMAN
Immunogen	Synthetic Peptide of Histone H2B (Tri Methyl Lys5)
Immunogen Region	
Specificity	The antibody detects endogenous Histone H2B (Tri Methyl Lys5) protein.
Immunogen Sequence	



Western blot analysis of 1) HeLa, 2) Rat Brain, 3) 3T3, diluted at 1:2000. Secondary antibody (STJ99666) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-503, Inventibiotec, MN, USA).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081