

Anti-Acetyl-Histone H3-Lys36 antibody (STJ197307)

STJ197307

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Histone H3-Lys36 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:200-500
Formulation	Liquid in PBS pH7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID [126961/333932/653604](#)
[3020/3021](#)
[H3C15.H3C14.H3C13](#)
[H3-3A.H3-3B](#)
[H32_HUMAN](#)
[H33_HUMAN](#)
[H31_HUMAN](#)

<

Immunogen	Synthetic Peptide of Histone H3 (Acetyl Lys36)
Immunogen Region	
Specificity	The antibody detects endogenous Histone H3 (Acetyl Lys36) protein.
Immunogen Sequence	



Western blot analysis of extracts from HeLa cells, untreated (-) or treated, 1:2000. Secondary antibody (STJ99666) was diluted at 1:20000 cells nucleus extracted by Minute TM Cytoplasmic and Nuclear Fractionation kit (SC-003, Inventiotech, MN, USA).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081