

Anti-Phospho-Histone H2A-Ser129 antibody (STJ197461)

STJ197461

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Histone H2A-Ser129 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS pH7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	723790/8337 8329/8330/8332/8336/8969 H2AC18.H2AC19 H2AC11.H2AC13.H2AC15.H2AC
Uniprot ID	H2A2A_HUMAN H2A1_HUMAN H2A3_HUMAN
Immunogen	Synthetic Peptide of Histone H2A (Phospho Ser129)
Immunogen Region	
Specificity	The antibody detects endogenous Histone H2A (Phospho Ser129) protein.
Immunogen Sequence	



Western blot analysis of extracts from HeLa cells, untreated (-) or treated (+) with Nocodazole, 1:5000. Secondary antibody (STJ99666) was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081