

Anti-SVIP antibody (26-76 aa) (STJ197512)

STJ197512

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Small VCP/p97-interacting protein (26-76 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	258010
Gene Symbol	SVIP
Uniprot ID	SVIP_HUMAN
Immunogen	Synthesized peptide derived from the human SVIP at the amino acid range 26-76
Immunogen Region	26-76 aa
Specificity	This antibody detects endogenous levels of human SVIP
Immunogen Sequence	