

Anti-SEMA3D antibody (442-492 aa) (STJ197542)

STJ197542

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Semaphorin-3D (442-492 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	223117
Gene Symbol	SEMA3D
Uniprot ID	SEM3D_HUMAN
Immunogen	Synthesized peptide derived from the human SEM3D at the amino acid range 442-492
Immunogen Region	442-492 aa
Specificity	This antibody detects endogenous levels of human SEM3D
Immunogen Sequence	