

**Anti-KIAA0825 antibody (199-249 aa) (STJ197548)**

STJ197548

**GENERAL INFORMATION**

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Uncharacterized protein KIAA0825 (199-249 aa) is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB   |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human/Mouse  |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 1 mg/mL   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:1000-2000  |
| <b>Formulation</b>         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

**TARGET INFORMATION**

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | <a href="#">285600</a>   |
| <b>Gene Symbol</b>        | <a href="#">KIAA0825</a>   |
| <b>Uniprot ID</b>         | <a href="#">K0825_HUMAN</a>  |
| <b>Immunogen</b>          | Synthesized peptide derived from the human K0825 at the amino acid range 199-249 |
| <b>Immunogen Region</b>   | 199-249 aa   |
| <b>Specificity</b>        | This antibody detects endogenous levels of human K0825                           |
| <b>Immunogen Sequence</b> |  |