

Anti-H2AC20 antibody (70-120 aa) (STJ197583)

STJ197583

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Histone H2A type 2-C (70-120 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8338
Gene Symbol	H2AC20
Uniprot ID	H2A2C_HUMAN
Immunogen	Synthesized peptide derived from the human H2A2C at the amino acid range 70-120
Immunogen Region	70-120 aa
Specificity	This antibody detects endogenous levels of human H2A2C
Immunogen Sequence	