

Anti-HS3ST1 antibody (201-251 aa) (STJ197590)

STJ197590

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Heparan sulfate glucosamine 3-O-sulfotransferase 1 (201-251 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9957
Gene Symbol	HS3ST1
Uniprot ID	HS3S1_HUMAN
Immunogen	Synthesized peptide derived from the human HS3S1 at the amino acid range 201-251
Immunogen Region	201-251 aa
Specificity	This antibody detects endogenous levels of human HS3S1
Immunogen Sequence	