

Anti-MTRNR2L3 antibody (1-24 aa) (STJ197606)

STJ197606

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Humanin-like 3 (1-24 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	MTRNR2L3
Uniprot ID	HMN3_HUMAN
Immunogen	Synthesized peptide derived from the human HMN3 at the amino acid range 1-24
Immunogen Region	1-24 aa
Specificity	This antibody detects endogenous levels of human HMN3
Immunogen Sequence	