

Anti-EXOC3L2 antibody (271-321 aa) (STJ197615)

STJ197615

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Exocyst complex component 3-like protein 2 (271-321 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	EXOC3L2
Uniprot ID	EX3L2_HUMAN
Immunogen	Synthesized peptide derived from the human EX3L2 at the amino acid range 271-321
Immunogen Region	271-321 aa
Specificity	This antibody detects endogenous levels of human EX3L2
Immunogen Sequence	