

Anti-HERC1 antibody (509-559 aa) (STJ197623)

STJ197623

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Probable E3 ubiquitin-protein ligase HERC1 (509-559 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8925
Gene Symbol	HERC1
Uniprot ID	HERC1_HUMAN
Immunogen	Synthesized peptide derived from the human HERC1 at the amino acid range 509-559
Immunogen Region	509-559 aa
Specificity	This antibody detects endogenous levels of human HERC1
Immunogen Sequence	