

Anti-C2orf76 antibody (29-79 aa) (STJ197629)

STJ197629

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-UPF0538 protein C2orf76 (29-79 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	130355
Gene Symbol	C2orf76
Uniprot ID	CB076_HUMAN
Immunogen	Synthesized peptide derived from the human CB076 at the amino acid range 29-79
Immunogen Region	29-79 aa
Specificity	This antibody detects endogenous levels of human CB076
Immunogen Sequence	