

Anti-AFF2 antibody (1104-1154 aa) (STJ197643)

STJ197643

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-AF4/FMR2 family member 2 (1104-1154 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2334
Gene Symbol	AFF2
Uniprot ID	AFF2_HUMAN
Immunogen	Synthesized peptide derived from the human AFF2 at the amino acid range 1104-1154
Immunogen Region	1104-1154 aa
Specificity	This antibody detects endogenous levels of human AFF2
Immunogen Sequence	