

Anti-LTBP2 antibody (1080-1130 aa) (STJ197655)

STJ197655

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Latent-transforming growth factor beta-binding protein 2 (1080-1130 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	4053
Gene Symbol	LTBP2
Uniprot ID	LTBP2_HUMAN
Immunogen	Synthesized peptide derived from the human LTBP2 at the amino acid range 1080-1130
Immunogen Region	1080-1130 aa
Specificity	This antibody detects endogenous levels of human LTBP2
Immunogen Sequence	