

Anti-TBC1D32 antibody (202-252 aa) (STJ197668)

STJ197668

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein broad-minded (202-252 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at-20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	221322
Gene Symbol	TBC1D32
Uniprot ID	BROML_HUMAN
Immunogen	Synthesized peptide derived from the human BROMI at the amino acid range 202-252
Immunogen Region	202-252 aa
Specificity	This antibody detects endogenous levels of human BROMI
Immunogen Sequence	