

Anti-ATP5F1EP2 antibody (1-51 aa) (STJ197674)

STJ197674

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-ATP synthase subunit epsilon-like protein, mitochondrial (1-51 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	ATP5F1EP2
Uniprot ID	AT5EL_HUMAN
Immunogen	Synthesized peptide derived from the human AT5EL at the amino acid range 1-51
Immunogen Region	1-51 aa
Specificity	This antibody detects endogenous levels of human AT5EL
Immunogen Sequence	