

Anti-PRAM1 antibody (310-360 aa) (STJ197690)

STJ197690

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-PML-RARA-regulated adapter molecule 1 (310-360 aa) is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range | WB 1:1000-2000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 84106 |
| Gene Symbol | PRAM1 |
| Uniprot ID | PRAM_HUMAN |
| Immunogen | Synthesized peptide derived from the human PRAM at the amino acid range 310-360 |
| Immunogen Region | 310-360 aa |
| Specificity | This antibody detects endogenous levels of human PRAM |
| Immunogen Sequence | |