

Anti-PRUNE1 antibody (48-98 aa) (STJ197699)

STJ197699

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Exopolyphosphatase PRUNE1 (48-98 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	58497
Gene Symbol	PRUNE1
Uniprot ID	PRUN1_HUMAN
Immunogen	Synthesized peptide derived from the human PRUNE at the amino acid range 48-98
Immunogen Region	48-98 aa
Specificity	This antibody detects endogenous levels of human PRUNE
Immunogen Sequence	