

Anti-UNC13C antibody (1996-2046 aa) (STJ197730)

STJ197730

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein unc-13 homolog C (1996-2046 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	440279
Gene Symbol	UNC13C
Uniprot ID	UN13C_HUMAN
Immunogen	Synthesized peptide derived from the human UN13C at the amino acid range 1996-2046
Immunogen Region	1996-2046 aa
Specificity	This antibody detects endogenous levels of human UN13C
Immunogen Sequence	