

Anti-C19orf25 antibody (14-64 aa) (STJ197733)

STJ197733

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-UPF0449 protein C19orf25 (14-64 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	148223
Gene Symbol	C19orf25
Uniprot ID	CS025_HUMAN
Immunogen	Synthesized peptide derived from the human CS025 at the amino acid range 14-64
Immunogen Region	14-64 aa
Specificity	This antibody detects endogenous levels of human CS025
Immunogen Sequence	