

Anti-TRAT1 antibody (56-106 aa) (STJ197838)

STJ197838

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-T-cell receptor-associated transmembrane adapter 1 (56-106 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50852
Gene Symbol	TRAT1
Uniprot ID	TRAT1_HUMAN
Immunogen	Synthesized peptide derived from the human TRAT1 at the amino acid range 56-106
Immunogen Region	56-106 aa
Specificity	This antibody detects endogenous levels of human TRAT1
Immunogen Sequence	