

Anti-LRRC36 antibody (380-430 aa) (STJ197858)

STJ197858

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Leucine-rich repeat-containing protein 36 (380-430 aa) is suitable for use in Western Blot research applications. |
| Applications | WB |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 1 mg/mL |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Dilution Range | WB 1:1000-2000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|--|
| Gene ID | 55282 |
| Gene Symbol | LRRC36 |
| Uniprot ID | LRRC36_HUMAN |
| Immunogen | Synthesized peptide derived from the human LRC36 at the amino acid range 380-430 |
| Immunogen Region | 380-430 aa |
| Specificity | This antibody detects endogenous levels of human LRC36 |
| Immunogen Sequence | |