

Anti-BMERB1 antibody (152-202 aa) (STJ197893)

STJ197893

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-bMERB domain-containing protein 1 (152-202 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	89927
Gene Symbol	BMERB1
Uniprot ID	MERB1_HUMAN
Immunogen	Synthesized peptide derived from the human CP045 at the amino acid range 152-202
Immunogen Region	152-202 aa
Specificity	This antibody detects endogenous levels of human CP045
Immunogen Sequence	