

Anti-ADA2 antibody (243-293 aa) (STJ197901)

STJ197901

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Adenosine deaminase 2 (243-293 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51816
Gene Symbol	ADA2
Uniprot ID	ADA2_HUMAN
Immunogen	Synthesized peptide derived from the human CECR1 at the amino acid range 243-293
Immunogen Region	243-293 aa
Specificity	This antibody detects endogenous levels of human CECR1
Immunogen Sequence	