

Anti-THAP9 antibody (238-288 aa) (STJ197910)

STJ197910

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-DNA transposase THAP9 (238-288 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79725
Gene Symbol	THAP9
Uniprot ID	THAP9_HUMAN
Immunogen	Synthesized peptide derived from the human THAP9 at the amino acid range 238-288
Immunogen Region	238-288 aa
Specificity	This antibody detects endogenous levels of human THAP9
Immunogen Sequence	