

Anti-MZT2A antibody (72-122 aa) (STJ197930)

STJ197930

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Mitotic-spindle organizing protein 2A (72-122 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	653784
Gene Symbol	MZT2A
Uniprot ID	MZT2A_HUMAN
Immunogen	Synthesized peptide derived from the human MZT2A at the amino acid range 72-122
Immunogen Region	72-122 aa
Specificity	This antibody detects endogenous levels of human MZT2A
Immunogen Sequence	