

**Anti-HEATR4 antibody (502-552 aa) (STJ197971)**

STJ197971

**GENERAL INFORMATION**

|                          |  |
|--------------------------|--|
| <b>Product Type</b>      | Primary antibodies   |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-HEAT repeat-containing protein 4 (502-552 aa) is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB   |
| <b>Host/Source</b>       | Rabbit   |
| <b>Reactivity</b>        | Human/Rat/Mouse  |

**PRODUCT PROPERTIES**

|                            |   |
|----------------------------|---|
| <b>Clonality</b>           | Polyclonal  |
| <b>Clone ID</b>            |   |
| <b>Concentration</b>       | 1 mg/mL   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:1000-2000  |
| <b>Formulation</b>         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.   |
| <b>Isotype</b>             | IgG   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

**TARGET INFORMATION**

|                           |  |
|---------------------------|--|
| <b>Gene ID</b>            | 399671   |
| <b>Gene Symbol</b>        | HEATR4   |
| <b>Uniprot ID</b>         | HEAT4_HUMAN  |
| <b>Immunogen</b>          | Synthesized peptide derived from the human HEAT4 at the amino acid range 502-552 |
| <b>Immunogen Region</b>   | 502-552 aa   |
| <b>Specificity</b>        | This antibody detects endogenous levels of human HEAT4                           |
| <b>Immunogen Sequence</b> |  |