

Anti-HOXD4 antibody (99-149 aa) (STJ197984)

STJ197984

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Homeobox protein Hox-D4 (99-149 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3233
Gene Symbol	HOXD4
Uniprot ID	HXD4_HUMAN
Immunogen	Synthesized peptide derived from the human HXD4 at the amino acid range 99-149
Immunogen Region	99-149 aa
Specificity	This antibody detects endogenous levels of human HXD4
Immunogen Sequence	