

Anti-RNF19B antibody (272-322 aa) (STJ197996)

STJ197996

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E3 ubiquitin-protein ligase RNF19B (272-322 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	127544
Gene Symbol	RNF19B
Uniprot ID	RN19B_HUMAN
Immunogen	Synthesized peptide derived from the human RNF19B at the amino acid range 272-322
Immunogen Region	272-322 aa
Specificity	This antibody detects endogenous levels of human RNF19B
Immunogen Sequence	