

Anti-GOLGA6A antibody (123-173 aa) (STJ198014)

STJ198014

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Golgin subfamily A member 6A (123-173 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	342096
Gene Symbol	GOLGA6A
Uniprot ID	GOG6A_HUMAN
Immunogen	Synthesized peptide derived from the human GOG6A at the amino acid range 123-173
Immunogen Region	123-173 aa
Specificity	This antibody detects endogenous levels of human GOG6A
Immunogen Sequence	