

Anti-ZMYM1 antibody (904-954 aa) (STJ198020)

STJ198020

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Zinc finger MYM-type protein 1 (904-954 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	79830
Gene Symbol	ZMYM1
Uniprot ID	ZMYM1_HUMAN
Immunogen	Synthesized peptide derived from the human ZMYM1 at the amino acid range 904-954
Immunogen Region	904-954 aa
Specificity	This antibody detects endogenous levels of human ZMYM1
Immunogen Sequence	