

Anti-FAM131C antibody (130-180 aa) (STJ198034)

STJ198034

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Protein FAM131C (130-180 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	348487
Gene Symbol	FAM131C
Uniprot ID	F131C_HUMAN
Immunogen	Synthesized peptide derived from the human F131C at the amino acid range 130-180
Immunogen Region	130-180 aa
Specificity	This antibody detects endogenous levels of human F131C
Immunogen Sequence	