

Anti-RASGRP4 antibody (487-537 aa) (STJ198071)

STJ198071

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RAS guanyl-releasing protein 4 (487-537 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	115727
Gene Symbol	RASGRP4
Uniprot ID	GRP4_HUMAN
Immunogen	Synthesized peptide derived from the human GRP4 at the amino acid range 487-537
Immunogen Region	487-537 aa
Specificity	This antibody detects endogenous levels of human GRP4
Immunogen Sequence	