

Anti-EML6 antibody (1135-1185 aa) (STJ198080)

STJ198080

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Echinoderm microtubule-associated protein-like 6 (1135-1185 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	400954
Gene Symbol	EML6
Uniprot ID	EMAL6_HUMAN
Immunogen	Synthesized peptide derived from the human EMAL6 at the amino acid range 1135-1185
Immunogen Region	1135-1185 aa
Specificity	This antibody detects endogenous levels of human EMAL6
Immunogen Sequence	