

Anti-CTAGE6 antibody (169-219 aa) (STJ198092)

STJ198092

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-cTAGE family member 6 (169-219 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	340307
Gene Symbol	CTAGE6
Uniprot ID	CTGE6_HUMAN
Immunogen	Synthesized peptide derived from the human CTGE6 at the amino acid range 169-219
Immunogen Region	169-219 aa
Specificity	This antibody detects endogenous levels of human CTGE6
Immunogen Sequence	