

Anti-SPAG11B antibody (38-88 aa) (STJ198145)

STJ198145

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Sperm-associated antigen 11B (38-88 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10407/653423
Gene Symbol	SPAG11B
Uniprot ID	SG11B_HUMAN
Immunogen	Synthesized peptide derived from the human SG11B at the amino acid range 38-88
Immunogen Region	38-88 aa
Specificity	This antibody detects endogenous levels of human SG11B
Immunogen Sequence	