

Anti-MARCHF11 antibody (198-248 aa) (STJ198153)

STJ198153

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-E3 ubiquitin-protein ligase MARCHF11 (198-248 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	441061
Gene Symbol	MARCHF11
Uniprot ID	MARHB_HUMAN
Immunogen	Synthesized peptide derived from the human MARHB at the amino acid range 198-248
Immunogen Region	198-248 aa
Specificity	This antibody detects endogenous levels of human MARHB
Immunogen Sequence	