

Anti-C8orf76 antibody (149-199 aa) (STJ198154)

STJ198154

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Uncharacterized protein C8orf76 (149-199 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84933
Gene Symbol	C8orf76
Uniprot ID	CH076_HUMAN
Immunogen	Synthesized peptide derived from the human CH076 at the amino acid range 149-199
Immunogen Region	149-199 aa
Specificity	This antibody detects endogenous levels of human CH076
Immunogen Sequence	