

Anti-CCDC144A antibody (1116-1166 aa) (STJ198188)

STJ198188

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Coiled-coil domain-containing protein 144A (1116-1166 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9720
Gene Symbol	CCDC144A
Uniprot ID	C144A_HUMAN
Immunogen	Synthesized peptide derived from the human C144A at the amino acid range 1116-1166
Immunogen Region	1116-1166 aa
Specificity	This antibody detects endogenous levels of human C144A
Immunogen Sequence	