

Anti-NDST3 antibody (132-182 aa) (STJ198196)

STJ198196

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Bifunctional heparan sulfate N-deacetylase/N-sulfotransferase 3-glucosamine N-sulfotransferase
Applications	NDST3 (132-182 aa) is suitable for use in Western Blot research applications.
Host/Source	WB
Reactivity	Rabbit
	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000-2000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9348
Gene Symbol	NDST3
Uniprot ID	NDST3_HUMAN
Immunogen	Synthesized peptide derived from the human NDST3 at the amino acid range 132-182
Region	132-182 aa
Specificity	This antibody detects endogenous levels of human NDST3
Immunogen Sequence	