

Anti-SPINK7 antibody (1-50 aa) (STJ198359)

STJ198359

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine protease inhibitor Kazal-type 7 (1-50 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84651
Gene Symbol	SPINK7
Uniprot ID	ISK7_HUMAN
Immunogen	Synthesized peptide derived from the human ISK7 at the amino acid range 1-50
Immunogen Region	1-50 aa
Specificity	This antibody detects endogenous levels of ISK7 at Human, Mouse, Rat
Immunogen Sequence	



Western Blot analysis of HeLa 293T using primary antibody at 1:1000 dilution 4A°C, overnight. Secondary antibody (STJS000791) was diluted at 1:10000 25A°C, 1.5 hours

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081