

Anti-COMMD6 antibody (1-50 aa) (STJ198367)

STJ198367

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-COMM domain-containing protein 6 (1-50 aa) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000
Formulation	Liquid in PBS containing 50% Glycerol and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	170622
Gene Symbol	COMMD6
Uniprot ID	COMD6_HUMAN
Immunogen	Synthesized peptide derived from the human COMMD6 at the amino acid range 1-50
Immunogen Region	1-50 aa
Specificity	This antibody detects endogenous levels of COMMD6 at Human, Mouse, Rat
Immunogen Sequence	